healthsmart®

Non-Contact Infrared Thermometer

Item #18-545-000

Please read this guidebook
completely before operating this unit.

Includes: 1 Thermometer, 1 User's Manual, 2 AAA Batteries

Index

Introduction	2
Intended for Use	3
Safety Precautions	3-4
Recommendations	5
Body Temperature Chart	5
Product Overview / LCD Display	6
Temperature Taking Hints	7
Before Use	7
Changing Measurement Modes (Forehead & Object)	8
Switching Between Fahrenheit and Celsius	8
How To Take A Reading - Forehead	9
How To Take A Reading - Object	10
Memory Recall	11
Replacing The Battery	12
FCC Statement and EMC Compatibility	12-13
Maintenance, Storage and Calibration	13
Troubleshooting	14 - 15
Product Specifications	16
Normalized Symbols	16-17
Important Information Regarding Electromagnetic Compatibility (EMC)	17 - 23
Warranty	23

Introduction

- Thank you for purchasing the Non-Contact Infrared Thermometer.
- Remove safety tab from battery. Open battery case and gently pull out the plastic safety tab from underneath the battery. Thermometer will not activate until this has been performed.
- For safe and proper use of this product, be sure to read and fully understand Safety Precautions contained in this instruction manual.
- Keep this Manual with you for a convenient reference.
- The reference body site is core and mode of operation is adjusted mode.
- The user can measure and change the battery under normal circumstances and maintain the device and its accessories according to the user manual.

Intended for Use

The device is intended for the intermittent measurement and monitoring of human body temperature by way of contact-free measurements (at a distance of approximately 0"~2" or 0-5 cm) on the forehead and object by consumers in the home. This device is intended for use on people of all ages.

Safety Precautions

Proper technique is critical to obtaining accurate temperatures. Therefore, in order to obtain an accurate measurement, please read this manual carefully before using.

- Use of this thermometer is not intended as a substitution for consultation with your physician. Measurement results are for reference only. Contact your physician if you have or suspect any health concerns.
- The infrared sensor must be kept clean, dry and undamaged at all times to ensure accurate measurements.
- Oils or lotions on the forehead may cause an inaccurate measurement, so make sure the forehead measurement site is clean.
- Do not touch or blow on the infrared sensor. A dirty sensor may cause inaccuracy.
- If the temperature of the storage area differs greatly from that of the measurement area, please wait for the thermometer temperature to equalize to the room temperature for about 30 minutes before use.
- Keep the device out of the reach of children and pets to avoid inhalation or swallowing of small parts. Do not allow children to take their own temperatures unattended. Children may not be able to use the device according to the instructions in this user manual. This is not a toy.
- Do not discard the battery into a fire. The battery is likely to explode.
- Do not store the thermometer in extreme temperatures below -4°F (-20°C) or over 131°F (55°C) or in humidity below 15% RH or over 90% RH. Failing to do so may cause inaccuracy.

- Do not use the thermometer if the main body is damaged (for example the infrared sensor is broken) The continuous use of a damaged unit may result in false readings or injury.
- Do not take the thermometer apart, repair or change any parts on the unit at any time.
- Do not use a mobile phone nearby when the thermometer is working.
- If your thermometer will not be used for a long time, it is advised to remove the battery to prevent possible damage due to chemical leakage. If the battery does leak, remove them carefully. Do not allow bare skin to touch the leaking fluid.
- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided. This unit can be affected by portable and mobile RF communications equipment.
- The materials (e.g. ABS) of expect contact with patient had passed the ISO 10993-5 and ISO 10993-10 standards test, no toxicity, allergy and irritations reaction. However, based on current science and technology, other potential allergic reactions are unknown. If you have allergic reaction to the materials, please stop using immediately and consult your physician.
- Store the thermometer in a dry, clean place. Do not leave the thermometer exposed to any chemical solvent, lint, dust direct sunshine or high temperatures.
- When using, do not touch the battery and the patient simultaneously.
- Do not maintain or service the device while the device is in use.

CAUTION: This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Recommendations

- Don't use the product for other purposes than what it is intended to be used.
- It is not a waterproof device. Clean it with a soft, wet cloth or cotton swab.
- Do not leave the product exposed to any chemical solvent, direct sunshine or high temperature.
- Don't drop or impose any vibration or impact on the product.
- Don't use a mobile phone nearby when the thermometer is working.
- Please act according to the native law in handling and disposing of the batteries.
- Take out the batteries if you are not going to use the unit for a long period of time.

Body Temperature Chart


A person's normal temperature will vary and trends to decrease with age. The following table shows normal temperature ranges by age:

Normal temperature ranges by age:


Age	Temperature (°F)	Temperature (°C)
0 - 2 years	97.5 - 100.4 °F	36.4 - 38.0 °C
3 - 10 years	97.0 - 100.0 °F	36.1 - 37.8 °C
11 - 65 years	96.6 - 99.7°F	35.9 - 37.6 °C
> 65 years	96.4 - 99.5°F	35.8 - 37.5 °C

The range of normal temperature varies from person to person and can be influenced by many factors such as time of day, measure from different sites, level of activity, medications, emotion, etc. We recommend that you practice with the thermometer when you or family members are healthy, so you are familiar with the readings when they are ill.

Product Overview

LCD Display

This package contains: (1) Thermometer, (2) 1.5v AAA Batteries, (1) User Manual

Temperature Taking Hints

1. External factors may influence forehead temperature, particularly when an individual has:
 - been exposed to very hot or very cold temperatures
 - recently been swimming or bathing
 - had their forehead coveredIn any of these cases listed above, remove the person from the situation or uncover their head and wait 30 minutes prior to taking a temperature.
2. With the object measurement, only the surface temperature of the object is measured. It can be different from its internal temperature, especially if the surface is exposed to direct sunlight or dryness.

Before Use

1. To achieve an accurate measurement, it is important to check the infrared sensor and ensure it is clean before use. To clean the sensor, gently wipe its surface with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 20 minutes drying time before taking temperatures.
2. Oils or lotion on the forehead may cause an inaccurate measurement, so make sure the forehead is clean. If you clean the forehead, wait 5~10 minutes before measuring.
3. If the temperature of the storage area differs greatly from that of the measuring area, please wait 30 minutes for the thermometer and room temperature to equalize before measuring.
4. Open the battery cover and remove the insulating plastic from the battery compartment. The thermometer should turn on automatically.

Changing Measurement Modes

NOTE: The Forehead measurement is the thermometer's default measuring mode.

To Measure in Forehead Mode:

1. Press the SCAN button to turn the thermometer ON. It will now be in Standby mode measuring ambient temperature.
2. Because the Forehead Mode is thermometer's default measuring mode, the thermometer is ready to take a forehead temperature measurement.
3. See "How To Take A Reading - Forehead" for more information.

To Measure in Object Mode:

1. Press the SCAN button to turn the thermometer ON. It will now be in Standby mode measuring ambient temperature.
2. Press the SCAN button and M (Memory) button together to switch from Standby mode to Object mode.
3. See "How To Take A Reading - Object" for more information.
4. Press the M (Memory) button again to return to Standby mode.

Switching Between Fahrenheit and Celsius

This thermometer is set to measure in Fahrenheit, but can be switched to measure in Celsius.

To switch from Fahrenheit to Celsius

1. Press the SCAN button to turn the thermometer ON. It will now be in Standby mode measuring ambient temperature.
2. Press and hold the M (Memory) button for 2 seconds until the screen changes from °F to °C. Once the thermometer beeps, the unit is ready for measurement.

NOTE: To switch the mode again, wait until the unit beeps and then press and hold the M (Memory) button again until the mode switches.

How To Take A Reading - Forehead

NOTE: Make sure the measurement sensor is clean before taking a reading. Be sure the thermometer is set to measure the temperature of a person. Attempting to take a measurement anywhere else on the body may result in an inaccurate reading.

1. If the thermometer is OFF, press any button to turn it ON. The screen should display the ambient temperature, *Fig 1*.


FIG. 1

2. Aim the front of the thermometer at the middle of the person's forehead, *Fig. 2*. Hold the thermometer approximately 0" ~ 2" (0~5 cm) away from the forehead.

FIG. 2

3. Press and hold the SCAN button. The screen should appear like *Fig. 3* and is ready to begin taking a measurement.

4. While continuing to hold the SCAN button, slowly move the thermometer toward or away from the person's forehead. You will hear a 'beeping' sound until the thermometer is at the correct distance. Once the beeping stops, hold the thermometer still as it is taking the measurement.

5. The thermometer will make a final 'beep', indicating the measurement is complete and display the result on the backlit LCD screen.


FIG. 3

6. After the test is complete, press the M (Memory) button or wait 1 minute for the display to return to standby mode and display the ambient temperature. After 1 minute in standby mode, the thermometer will automatically shut OFF.

NOTE: The LED light around the SCAN button will illuminate according to the temperature reading.

- GREEN = 93.2°F ~ 100.3°F (34.0°C ~ 37.9°C)
- RED = 100.4°F ~ 109.3°F (38.0°C ~ 42.9°C)
- LO = < 89.6°F (32.0°C)
- HI = > 109.3°F (42.9°C)

How To Take A Reading - Object

NOTE: Make sure the measurement sensor is clean before taking a reading. Be sure the thermometer is set to measure the temperature of an object.

1. If the thermometer is off, press any button to turn it ON. The screen should display the ambient temperature, *Fig 1*.

FIG. 1

2. Press the SCAN button and M (Memory) button together to switch between Object and Forehead mode. The screen should appear like *Fig. 2* when in Object mode.

FIG. 2

3. Aim the front of the thermometer at the object you want to measure, *Fig. 3*. Hold the thermometer approximately 0" ~ 2" (0~5 cm) away from the object.

FIG. 3

4. Press and hold the SCAN button. The screen should appear like *Fig. 4* and begin taking the measurement.

NOTE: Holding the measurement button down will give a continuous reading until the button is released.

5. After the test is complete, press the M (Memory) button or wait 1 minute for the display to return to standby mode and display the ambient temperature. After 1 minute in standby mode, the thermometer will automatically shut OFF.

NOTE: The LED light will NOT illuminate as it does in the Forehead measurement mode.

FIG. 4

Memory Recall

There are a total of 9 memory readings.

NOTE: Only Forehead measurements are recorded. Object and Ambient temperature measurements are NOT stored in the memory.

The current measurement is always stored in the last position (number 9).

To recall measurements in memory:

1. If the thermometer is OFF, press any button to turn it ON.
2. Press the M (Memory) button. It will show the position of the most current measurement *Fig 1*. Followed by the actual temperature measurement *Fig 2*.

NOTE: All subsequent measurements taken will fill the memory spaces until they reach the number 1 spot. Once all of the spaces are filled, the oldest measurement will automatically be deleted from memory.

3. Continue to press the M (Memory) button to scroll the stored measurements.
4. Press the SCAN button and M (Memory) button at the same time to return to the Standby mode.

FIG. 1

FIG. 2

Replacing the Battery

When the low battery indicator “” appears on the screen, or the thermometer does not turn ON, it is time to replace the batteries.

To Replace the batteries:

1. Gently slide the battery cover back.
2. Carefully remove the old batteries and properly discard.
3. Insert 2 new AAA 1.5v batteries. Ensure they are inserted correctly according to the proper polarity,
Fig 1.


FIG. 1

4. Slide the battery cover back on until it snaps into place.

To protect the environment, please dispose of exhausted batteries according to current federal, state and local regulations. Keep batteries out of the reach of children.

FCC Statement and EMC Compatibility

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Statement and EMC Compatibility

15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Maintenance, Storage and Calibration

Maintenance:

1. To ensure accurate measurements, keep the infrared sensor and distance sensor tip clean and free of scratches. Finger prints or dirt will affect the accuracy of the thermometer.
2. To clean the infrared sensor and distance sensor: gently wipe the sensors with a cotton swab slightly moistened with alcohol and immediately wipe dry with a clean cotton swab (Fig 1). The thermometer can be used after the alcohol dries.
3. Use a soft dry cloth to clean the body of the thermometer (Fig 2).
4. Never clean the thermometer with an abrasive cleanser, thinner, benzene or submerge the thermometer into water or other liquids.

FIG. 1

FIG. 2



Storage:

1. After each use clean the thermometer.
2. Store in a clean dry place at room temperature. NEVER leave the thermometer exposed to direct sun, high temperatures, humidity, fire, flames, shocks or hits.
3. If you do not use the thermometer for a long time, take out the batteries.

Calibration:

1. The thermometer is factory-calibrated. If the thermometer is used according to the user manual, periodic calibration is not required. Do not try to modify or reassemble the thermometer.

Troubleshooting

Problem	Check lists	Action
No response/ Automatic reset when the plastic battery insulator tab is pulled out	Battery used up?	Change new battery
	Battery in wrong polarity?	Take out battery, reinsert battery correctly
	Poor battery contact	
	The measured temperatures is lower than 89.6°F (32.0°C) or higher than 109.3°F (42.9°C). Check the method of operation	Follow the instructions for use.
	Indicates a hardware problem	Contact your distributor
	Operating temperatures is out of range	Use the thermometer within the range of operating temperatures
The measurement is not accurate or if there is any doubt on the measured result	Please check if the infrared sensor is clean	Clean the infrared sensor with a cotton swab according to the instructions.
	Please check if the measurement method is correct and in the correct mode	Ensure you have fully read the user instructions and know how to properly use the thermometer.

Troubleshooting

Problem	Check lists	Action
The measurement is not accurate or if there is any doubt on the measured result	Please check if you have let the thermometer and patient temperature get stabilized in the room for at least 30 minutes	Please keep the thermometer and patient in the measuring room at least 30 minutes before using
	Are you using the thermometer indoor?	Please take the measurement indoor
	Please check if you held the thermometer in your hand too long and affect the accuracy	Put the thermometer on the table in the room where the measurement is taking place and let it cool down first
	Please check if there is low battery symbol on LCD	Change new battery

Specifications

- **Product Name:** Non-Contact Infrared Thermometer
- **Power supply:** DC 3V (2 x AAA battery)
- **Forehead Measuring Range:** 89.6°F – 109.3°F (32.0°C – 42.9°C)
Forehead Measuring Accuracy: $\pm 0.4^{\circ}\text{F}$ (or $\pm 0.2^{\circ}\text{C}$)
from 95°F – 107.6°F (or 35°C – 42.0°C) $\pm 0.5^{\circ}\text{F}$ ($\pm 0.3^{\circ}\text{C}$) for other temperature range
- **Object Measuring Range:** -7.6°F – 176.0°F (-22°C – 80°C)
- **Object Measuring Accuracy:** $\pm 4^{\circ}\text{F}$ ($\pm 2^{\circ}\text{C}$)
- **Ambient Temperature Range:** 41.0°F – 139.8°F (5.0°C – 59.9°C)
- **Ambient Temperature Accuracy:** $\pm 2^{\circ}\text{F}$ ($\pm 1^{\circ}\text{C}$)
- **Measurement Time:** 1 Second
- **Measurement Interval:** 1 Second
- **Dimension:** 5.12 (L) x 1.77 (W) x 1.98(D) in. [130.1(L) x 45(W) x 50.3(D) mm]
- **Weight:** 72g (without batteries)
- **Operation Temperature:** 59.0°F – 104.0°F (15°C – 40°C)
- **Operation Relative Humidity:** 15% RH – 85%
- **Storage & Transportation Temperature:** -4°F – 131°F (-20°C – 55°C)
- **Storage & Transportation Relative Humidity:** 15% RH – 90% RH
- **Service Life:** 2 years
- **Battery Life:** Approximately 6 months with 5 tests per day

Normalized Symbols

Applied part of type BF

Disposal in accordance with Directive
2012/19/EU (WEEE)

The name and the address of the manufacturer.

Caution

Refer to Instruction Manual.

Normalized Symbols

IP22 The first number 2: Protected against solid foreign objects of 12,5 mm Φ and greater. The second number: Protected against vertically falling water drops when enclosure tilted up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is tilted at any angle up to 15° on either side of the vertical.

Important Information Regarding Electromagnetic Compatibility (EMC)

With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation. Medical devices should also not interfere with other devices.

In order to regulate the requirements for EMC (Electro Magnetic Compatibility) with the aim to prevent unsafe product situations, the EN60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices. Medical devices manufactured by HealthSmart International conform to this EN60601-1-2:2007 standard for both immunity and emissions.

Nevertheless, special precautions need to be observed: The use of accessories and cables other than those specified by manufacturer, with the exception of cables sold by manufacturer as replacement parts for internal components, may result in increased emission or decreased immunity of the device.

- The medical devices should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the medical device should be observed to verify normal operation in the configuration in which it will be used.
- Refer to further guidance below regarding the EMC environment in which the device should be used.

Guidance and Manufacturer's Declaration

(Table 1)

Guidance and manufacturer's declaration – electromagnetic emissions		
The thermometer are intended for use in the electromagnetic environment specified below. The customer or the user of these thermometer should assure that it is used in such environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The thermometer use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The thermometer are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Not applicable	

(Table 2)

Guidance and manufacturer's declaration – electromagnetic immunity			
The thermometer are intended for use in the electromagnetic environment specified below. The customer or the user of these thermometer should assure that it is used in such environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floor should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply IEC 61000-4-11	Not applicable	Not applicable	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

(Table 4)

Guidance and manufacturer's declaration – electromagnetic immunity			
The thermometer are intended for use in the electromagnetic environment specified below. The customers or the users of these thermometer should assure that it is used in such environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	Not applicable		Portable and mobile RF communications equipment should be used no closer to any part of the thermometer including cables, than the recommended separation distance calculated from the equation appropriate to the frequency of the transmitter. Recommend separation distance $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey,

			<p>a. should be less than the compliance level in each frequency range.</p> <p>b. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	--	--	---

Note1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note2: These guidelines may not apply in all situations.

Electromagnetic propagation is

affected by absorption and reflection from structures, objects, and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the thermometer are used exceeds the applicable RF compliance level above, the thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the thermometer.

^b Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

(Table 6)

Recommended separation distance between portable and mobile RF communications equipment and the thermometer			
<p>The thermometer are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customers or the users of these thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the thermometer as recommended below, according to the maximum output power of the communications equipment.</p>			
Output Power of Transmitter in Watt	Separation distance according to frequency of transmitter in meter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>Note: At 80MHz and 800MHz, the separation distance for the higher frequency range applies</p> <p>Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</p>			

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Warranty

Three-year warranty is available from purchasing date, excluding user-caused failures listed below:

- (1) Failure resulted in unauthorized disassembly and modification
- (2) Failure resulted in unexpected drop during application or transportation
- (3) Failure resulted in operation away from proper instruction in user's manual

Declaration of Conformity

IEC60601-1, IEC60601-1-2, IEC60601-1-4, ASTM E1965-2009 revised, ISO10993-5, ISO10993-10, ISO10993-1, ISO 14971, EN 12470-5.

NOTES:

- 1) Under the environment with electrostatic discharge, the unit may malfunction and may require user to reset the unit.
- 2) Dispose with care, please consult the retailer for details.

hs
**Satisfaction
Guaranteed**

Manufactured for:

HealthSmart® International

1931 Norman Drive
Waukegan, IL 60085

Register this product at LiveHealthSmart.com

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©2016 HealthSmart® International