

# INSTRUCTIONS



## Bench-Top Fume Hood

Catalog No. 50000-0002

The molded polyethylene bench-top fume hood has two stack diameters, 4" (102mm) and 6" (152mm). If the 6" diameter stack is to be used, it is necessary to cut off the 4" diameter extension to avoid restricting the airflow. A small hand saw or sabre saw can be used for this purpose.



The sliding, clear-acrylic door has two black knobs (A) to raise and lower it. The side-tracks have three holes that are provided for the detent pins (B) which hold the door in position at the selected level.

When the detent pins are in the lowest position, a 1/4" wide gap will be left beneath the door. This gap will allow air to enter and sweep the deck. The detent pins should not be removed completely but should always be in one set of holes.

An exhaust system capable of drawing 200 cfm (free air) is recommended for use with this fume hood.

Do not place very hot items directly in contact with the plastic material. A practical temperature limit is 80°C (176°F). Do not use abrasive cleaners on the clear door as it will scratch and become cloudy.

If you require further assistance in the installation, operation, or use of this Scienceware® product, please contact us.

### BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA  
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

95000-0011


Scienceware® is a reg. TM of Bel-Art Products • Draft Checker™ is a trademark of Bel-Art Products  
Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.

03/02

## Draft Checker™ Airflow Detection Powder

## Safely Test Fume Hood Exhaust

Catalog No. 19939-0000



- Amorphous Silica Powder is inert, Non-Toxic, and Non-Flammable.
- Great for Detecting Drafts, Gas Leaks and Air Currents in Fume Hoods.
- 1.5gm Bottle.

**Just a squeeze or two from the plastic bottle dispenses a cloud of ultra fine lightweight silica powder.**

### BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA  
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

Contact your favorite dealer or visit our web site @ [belart.com](http://belart.com) to see our complete catalog and to request your own FREE copy of our Scienceware® catalog or call 1-800-4BEL-ART.